



TRANSFER SETTINGS

MiniMed™ 530G pump to MiniMed™ 770G pump

STEP 1: Copy your current settings onto the Insulin Pump Settings Form. If for any reason you cannot copy the settings, please contact your healthcare professional.

STEP 2: Insert a AA battery into your new pump. You can use either a lithium, alkaline, or rechargeable battery.

STEP 3: Follow the instructions found in this guide to enter the settings that you have recorded on the Insulin Pump Settings Form into your new pump. For additional information on your pump, please consult the MiniMed™ 530G System User Guide.

Medtronic

We recommend that you record all settings on this form to refer to as you enter them into the new pump.

BASAL SETTINGS

STANDARD BASAL RATE:

Main ▶ Basal ▶ Basal Review

- 1: 12:00AM _____ U/H
- 2: _____ U/H
- 3: _____ U/H
- 4: _____ U/H
- 5: _____ U/H
- 6: _____ U/H

MAX BASAL RATE:

Main ▶ Basal ▶ Basal Setup ▶ Max Basal Rate

_____ U/H

BASAL PATTERNS:

Main ▶ Basal ▶ Basal Setup ▶ Patterns

On Off

Basal Pattern A:

Main ▶ Basal ▶ Basal Setup ▶ Set/Edit Basal ▶ Pattern A

- 1: 12:00AM _____ U/H
- 2: _____ U/H
- 3: _____ U/H
- 4: _____ U/H

Basal Pattern B:

Main ▶ Basal ▶ Basal Setup ▶ Set/Edit Basal ▶ Pattern A

- 1: 12:00AM _____ U/H
- 2: _____ U/H
- 3: _____ U/H
- 4: _____ U/H

UTILITIES

ALARM

Main ▶ Utilities ▶ Alarm ▶

Alert Type: Beep Long Beep Medium
 Beep Short Vibrate

Auto Off: On Off _____ Hrs

Low Resv Warning: _____ Units _____ Time

BOLUS SETTINGS

BOLUS WIZARD SETUP:

Main ▶ Bolus ▶ Bolus Setup ▶ Bolus Wizard Setup ▶ Review Settings

Wizard: On Off **Carb Units:** Grams Exch

Carb Ratios: 1: 12:00AM _____ grams/U
 2: _____ grams/U
 3: _____ grams/U

Sensitivity: 1: 12:00AM _____ mg/dL/U
 2: _____ mg/dL/U
 3: _____ mg/dL/U

BG Target: 1: 12:00AM _____ - _____ mg/dL
 2: _____ - _____ mg/dL
 3: _____ - _____ mg/dL

Active Ins Time: _____ Hrs

MAX BOLUS:

Main ▶ Bolus ▶ Bolus Setup ▶ Max Bolus

_____ units

SCROLL RATE:

Main ▶ Bolus ▶ Bolus Setup ▶ Scroll Rate

0.1 U 0.05 U 0.025 U

DUAL/SQUARE BOLUS:

Main ▶ Bolus ▶ Bolus Setup ▶ Dual Square Bolus

On Off

EASY BOLUS:

Main ▶ Bolus ▶ Bolus Setup ▶ Easy Bolus

On/Set Off _____ units

BG REMINDER:

Main ▶ Bolus ▶ Bolus Setup ▶ BG Reminder

On Off

MISSED BOLUS REMINDER:

Main ▶ Bolus ▶ Bolus Setup ▶ Missed Bolus Reminder

On Off

1: _____
 2: _____

CONTINUOUS GLUCOSE MONITORING INITIATION SETTINGS

THRESHOLD SUSPEND ALARM (MINIMED® 530G SYSTEM ONLY)

On _____ mg/dL (60–90 mg/dL) Off

LOW GLUCOSE ALERT

Start on _____ (date)		
	Time	Units
	12 a.m.	_____ mg/dL
<input type="checkbox"/> On	_____	_____ mg/dL
	_____	_____ mg/dL
OR		
		<input type="checkbox"/> Off

LOW PREDICTIVE ALERT

On _____ minutes (5–30 min) Off

FALL RATE ALERT

On _____ mg/dL/minute (1.1–5.0) Off

LOW REPEAT

_____ minutes (5 min–1 hour)

Low Repeat is turned on to the default of 20 minutes when any Low Alert setting is programmed.

HIGH GLUCOSE ALERT

Start on _____ (date)		
	Time	Units
	12 a.m.	_____ mg/dL
<input type="checkbox"/> On	_____	_____ mg/dL
	_____	_____ mg/dL
OR		
		<input type="checkbox"/> Off

HIGH PREDICTIVE ALERT

On _____ minutes (5–30 min) Off

RISE RATE ALERT

On _____ mg/dL/minute (1.1–5.0) Off

HIGH REPEAT

_____ minutes (5 min–3 hours)

Default is set at 1 hour. May need to increase to 2 hours to minimize alerts.

STARTUP WIZARD

Once you've inserted the battery, the pump powers up and the **Startup Wizard** begins, it will guide you through setting language, time format, and current time and date.

1. Select **English**.
2. When the **Select Time Format** screen appears, select a **12 Hour** or a **24 Hour** time format.
3. When **Enter Time** screen appears, the hour will be flashing. Press **^** or **v** to the correct hour and press **O**.
4. The minutes will be flashing. Press **^** or **v** to the correct minutes and press **O**.
5. The **AM / PM** will be flashing. Press **^** or **v** if needed and press **O**. Select **Next**.
6. When the **Enter Date** screen appears, select **Year**. Press **^** or **v** to the correct year and press **O**. Repeat this step to set up **Month** and **Day**. Select **Next**.
7. The rewinding screen will briefly appear, followed by a confirmation message. Select **Ok**.



SETTING MAX BASAL

Menu > Options > Delivery Settings

1. Select **Max Basal/Bolus** and press **O**.
2. Select **Max Basal** to enter the **U/hr** amount. Use **^** and **v** to set the **Max Basal** amount. Press **O**.
3. Select **Save**.



SETTING MAX BOLUS

Menu > Options > Delivery Settings

1. Select **Max Basal/Bolus** and press **O**.
2. Select **Max Bolus** to enter the **U/hr** amount. Use **^** and **v** to set the **Max Bolus** amount. Press **O**.
3. Select **Save**.



SETTING BASAL PATTERN

Menu > Basal > Delivery Settings > Basal Pattern Setup

1. Select **Basal 1**.
2. Select **Options**, then select **Edit**.
3. Press **O** to enter **End** time. Press **O** to leave the **End** time at **12:00A**. If you are setting multiple rates use **^** or **v** to change the **End** time and press **O**.
4. Press **^** to enter **U/hr** and press **O**.
5. Continue adding **End** times and basal rates as needed, otherwise select **Done**.
6. Review the basal rates. If you need to correct, press **⬅** to go back to edit mode. If Basal rates are correct select **Save**.



TURN BOLUS WIZARD ON AND SETUP

Menu > Options > Delivery Settings > Bolus Estimate Setup

1. Select **Bolus Wizard** to turn **On**.
2. Press **✓** to see the rest of the message then select **Next**.





3. Review the description of **Carb Ratio** and select **Next**.
4. Press **○** to enter **End** time. Press **○** again to leave the **End** time at **12:00A**. If you are setting multiple Carb Ratios use **▲** or **▼** to change the **End** time, then press **○**.
5. Press **▲** or **▼** to enter g/U then press **○**.
Continue adding **End** times and **Carb Ratios** as needed.
6. Select **Next** when done.

7. Review the description of **Sensitivity Factor** and select **Next**.
8. Press **○** to enter **End** time. Press **○** again to leave the **End** time at **12:00A**. If you are setting multiple **Sensitivity Factors**, use **▲** or **▼** to change the **End** time, then press **○**.
9. Press **▲** or **▼** to enter **mg/dL per U**, then press **○**.
Continue adding **End** times and **Sensitivity Factors** as needed.
10. Select **Next** when done.

11. Review the description of **BG Target** and select **Next**.
12. Press **○** to enter **End** time. Press **○** again to leave the **End** time at **12:00A**. If you are setting multiple **BG Targets**, use **▲** or **▼** to change the **End** time, then press **○**.
13. Press **▲** or **▼** to enter **Lo** and then press **○**.
14. Press **▲** or **▼** to enter **Hi** and then press **○**.
Continue adding **End** times and **BG Targets** as needed.
15. Select **Next** when done.

16. Review the description of **Active Insulin Time** and select **Next**.
17. Press **○** to enter **Duration** time.
18. Press **▲** or **▼** to enter **hr** then press **○**.
19. Select **Save**.

The **Bolus Wizard** setup is now complete.




Note: The settings shown are for illustration purposes only — your settings may be different.

TURN SENSOR ON

Menu > Options > Utilities > Sensor Settings

1. Select **Sensor** to turn feature **On**.
You can now access the **SmartGuard** menu and enter the settings.



PROGRAMMING HIGH SETTINGS

Menu > Options > SmartGuard

1. Select **High Setup**.
2. Press **○** on the time segment. Press **○** again to leave the **End** time at **12:00A**. If setting multiple time segments, press **▲** to enter **End** time of first segment, and press **○**.
3. Press **▲** or **▼** to set **Hi** limit and press **○**. Press **○** to continue.
4. Select each feature you wish to turn on. If a feature is on, select it again to turn it back off.
5. Once settings are selected, select **Next** and then select **Done**.
6. Verify that settings are correct and select **Save**.
7. If snooze time needs to be changed, press **▼** to **Snooze** and press **○**.
Select **High Snooze**. Press **▲** or **▼** to enter the correct time and press **○**. Select **Save**.




Your **High Setup** is now complete.

PROGRAMMING LOW SETTINGS

Menu > Options > SmartGuard

1. Select **Low Setup**.
2. Press **○** on the time segment. Press **○** again to leave the **End** time at **12:00A**. If setting multiple time segments, press **▲** to enter **End** time of first segment, and press **○**.
3. Press **▲** or **▼** to set **Lo** limit and press **○**. Press **○** to continue.
4. Select each feature you wish to turn on. If a feature is on, select it again to turn it back off.
5. Once settings are selected, select **Next** and then select **Done**.
6. Verify that settings are correct and select **Save**.
7. If snooze time needs to be changed, press **▼** to **Snooze** and press **○**.
Select **Low Snooze**. Press **▲** or **▼** to the correct time and press **○**. Select **Save**.



Your **Low Setup** is now complete.


Turn to page 5 to continue programming sensor settings.

TURN AUTO MODE ON

Menu > Options > SmartGuard

1. Select **Auto Mode**.
2. Select **Auto Mode** again to turn Auto Mode **On**.
3. Check the screen to make sure that Auto Mode is set to **On**.
Select **Save**

Note: If SmartGuard Auto Mode is not yet ready, after selecting **Save**, an alert occurs instructing to check the Auto Mode Readiness screen.




PAIR THE TRANSMITTER

Menu > Options > Utilities > Device Options > Pair Device

1. Attach your transmitter to the charger and place next to the pump.
2. Select **Search** on the pump and immediately remove the transmitter from the charger.
3. The **Select Device** screen appears with a list of available devices. Select the **CGM** device that matches the serial number (**SN**) on your transmitter.
4. Make sure the serial number on the pump screen and transmitter match and select **Confirm**.

Your pump displays a message if the pump and transmitter are paired successfully. If the Sensor feature is turned on, the Connection icon  appears on the Home screen.



Please take some time to review ALL your entries and ensure they are correct.

In Need of Technical Assistance?

24-Hour Technical Support Is Here to Help

Call **1.800.646.4633**

For medical advice, contact your healthcare team.

Note: If you turn Auto Mode on before the 48 hour warm up is complete, you may receive an alert at midnight to Enter a BG for Auto Mode. To avoid this, you may want to wait to turn on Auto Mode on day 3 of your pump wear.