



# TRANSFER SETTINGS

MiniMed™ 630G Pump to MiniMed™ 780G Pump

## Transfer Settings: MiniMed 630G pump to MiniMed 780G pump

**STEP 1** – Copy your current settings onto the insulin pump settings form. If for any reason you cannot copy the settings, please contact your healthcare professional.

**STEP 2** – Insert a AA battery into your new pump. You can use either a lithium, alkaline, or rechargeable battery.

**STEP 3** – Follow the instructions found in this guide to enter the settings that you have recorded on the insulin pump settings form into your new pump. For additional information on your pump, please consult the MiniMed pump system user guide.

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

## Settings from the MiniMed 630G insulin pump

We recommend that you record all settings on this form and use this form as reference when you enter your settings into the new pump

### Insulin Settings, Utilities, Sensor Settings

#### INSULIN SETTINGS

##### BOLUS WIZARD SETUP:

*Menu > Insulin Settings > Bolus Wizard Setup*

Bolus Wizard:  On  Off

##### Carb Ratio:

12:00 A \_\_\_\_\_ g/U  
 \_\_\_\_\_ g/U  
 \_\_\_\_\_ g/U

##### Insulin Sensitivity Factor:

12:00 A \_\_\_\_\_ mg/dL/U  
 \_\_\_\_\_ mg/dL/U  
 \_\_\_\_\_ mg/dL/U

##### BG Target:

12:00 A \_\_\_\_\_ - \_\_\_\_\_ mg/dL  
 \_\_\_\_\_ - \_\_\_\_\_ mg/dL  
 \_\_\_\_\_ - \_\_\_\_\_ mg/dL

Active Insulin Time: \_\_\_\_\_ hr

##### BASAL PATTERN SETUP:

*Menu > Insulin Settings > Basal Pattern Setup*

##### Basal 1:

12:00 A \_\_\_\_\_ U/hr  
 \_\_\_\_\_ U/hr  
 \_\_\_\_\_ U/hr

##### Basal Pattern Name:

12:00 A \_\_\_\_\_ U/hr  
 \_\_\_\_\_ U/hr  
 \_\_\_\_\_ U/hr

##### PRESET TEMP SETUP:

*Menu > Insulin Settings > Preset Temp Setup*

Preset Temp Name: \_\_\_\_\_

Type:  Rate \_\_\_\_\_ U/hr  Percent \_\_\_\_\_ %

Duration \_\_\_\_\_ hr

##### DUAL/SQUARE WAVE:

*Menu > Insulin Settings > Dual/Square Wave*

Dual:  On  Off

Square:  On  Off

#### INSULIN SETTINGS (CONT.)

##### PRESET BOLUS SETUP:

*Menu > Insulin Settings > Preset Bolus Setup*

Preset Bolus Name: \_\_\_\_\_

Bolus: \_\_\_\_\_ U Type: \_\_\_\_\_

\_\_\_\_\_ Now \_\_\_\_\_ Square

Duration \_\_\_\_\_ hr

##### BOLUS INCREMENT:

*Menu > Insulin Settings > Bolus Increment*

0.025 U  0.05 U  0.1 U

##### MAX BASAL/BOLUS:

*Menu > Insulin Settings > Max Basal/Bolus*

Max Basal: \_\_\_\_\_ U/hr Max Bolus: \_\_\_\_\_ U

##### EASY BOLUS:

*Menu > Insulin Settings > Easy Bolus*

On  Off Step Size \_\_\_\_\_ U

##### AUTO SUSPEND:

*Menu > Insulin Settings > Auto Suspend*

Alarm:  On  Off

Time: \_\_\_\_\_ hr

##### BOLUS SPEED:

*Menu > Insulin Settings > Bolus Speed*

Standard  Quick

#### UTILITIES

##### DISPLAY OPTIONS:

*Menu > Utilities > Display Options*

Brightness \_\_\_\_\_  Auto

Backlight \_\_\_\_\_ sec / min

##### REMOTE BOLUS:

*Menu > Utilities > Remote Bolus*

On  Off

##### CARB UNIT:

*Menu > Utilities > Carb Unit*

Grams  Exchanges

#### SENSOR SETTINGS

*Menu > Sensor Settings*

##### Sensor:

On  Off

### Sensor Settings

HIGH SETTINGS <i>Menu &gt; Sensor Settings &gt; High Settings</i>				
Time Period	High Limit	Alert before High	Alert on High	Rise Alert <sup>a</sup>
12:00 A - _____	_____ mg/dL	<input type="checkbox"/> On _____ min	<input type="checkbox"/> On	<input type="checkbox"/> On _____ mg/dL/min or ↑   ↑↑   ↑↑↑
_____ - _____	_____ mg/dL	<input type="checkbox"/> On _____ min	<input type="checkbox"/> On	<input type="checkbox"/> On _____ mg/dL/min or ↑   ↑↑   ↑↑↑
_____ - _____	_____ mg/dL	<input type="checkbox"/> On _____ min	<input type="checkbox"/> On	<input type="checkbox"/> On _____ mg/dL/min or ↑   ↑↑   ↑↑↑
High Snooze: _____ minutes				

<sup>a</sup> Option only available when the Rise Alert feature is set to on

LOW SETTINGS <i>Menu &gt; Sensor Settings &gt; Low Settings</i>				
Time Period	Low Limit	Alert before Low	Suspend on Low	Alert on Low
12:00 A - _____	_____ mg/dL	<input type="checkbox"/> On	<input type="checkbox"/> On	<input type="checkbox"/> On
_____ - _____	_____ mg/dL	<input type="checkbox"/> On	<input type="checkbox"/> On	<input type="checkbox"/> On
_____ - _____	_____ mg/dL	<input type="checkbox"/> On	<input type="checkbox"/> On	<input type="checkbox"/> On
Low Snooze: _____ minutes				


## Programming the MiniMed 780G insulin pump

### Startup Wizard

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

Once you've inserted the battery, the pump powers up and the **Startup Wizard** begins; it will guide you through setting language, time format, and current time and date.

1. Select **English**.
2. When the **Select Time Format** screen appears, select a **12 Hour** or a **24 Hour** time format.
3. When the **Enter Time** screen appears, the hour will be flashing. Press **^** or **v** to the correct hour and press **O**.
4. The minutes will be flashing. Press **^** or **v** to the correct minutes and press **O**.
5. The **AM / PM** will be flashing. Press **^** or **v** if needed and press **O**. Select **Next**.
6. When the **Enter Date** screen appears, select **Year**. Press **^** or **v** to the correct year and press **O**. Repeat this step to set up **Month** and **Day**. Select **Next**.
7. The **Rewinding** screen will briefly appear, followed by a confirmation message. Select **OK**.





### Setting Max Basal

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

#### Menu > > Delivery Settings

1. Select **Max Basal/Bolus** and press **O**.
2. Select **Max Basal**. A confirmation message appears. Select **Continue**.
3. Select **Max Basal** to enter the **U/hr** amount. Use **^** and **v** to set the **Max Basal** amount. Press **O**.
4. Select **Save**.





### Setting Max Bolus

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

#### Menu > > Delivery Settings

1. Select **Max Basal/Bolus** and press **O**.
2. Select **Max Bolus**. A confirmation message appears. Select **Continue**.
3. Select **Max Bolus** to enter the units (**U**) amount. Use **^** and **v** to set the **Max Bolus** amount. Press **O**.
4. Select **Save**.



### Setting a Basal Pattern

**WARNING:** The pump is intended to be used with a basal pattern. The basal pattern must be manually entered and saved into the pump. The pump will operate with a basal rate of 0.0 U/hr until a basal pattern is entered and saved. There is no reminder message to program basal rates. Consult a healthcare professional to determine what basal pattern is needed. For more information about basal patterns, consult the MiniMed 780G System User Guide.

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

Menu > ⚙ > Delivery Settings > Basal Pattern Setup



1. Select **Basal 1**.
2. Select **Options**, then select **Edit**.
3. Press **⊙** to enter **End** time. Press **⊙** to leave the end time at **12:00A**. If you are setting multiple rates use **▲** or **▼** to change the end time and press **⊙**.
4. Press **▲** to enter **U/hr** and press **⊙**.
5. Continue adding end times and basal rates as needed. The end time for your last rate must be **12:00A**.
6. The **Review** option appears only when the last end time in your basal pattern is set to **12:00A**. Select **Review** after all rates are entered to go to the next screen.

**Note:** Your basal rates have not yet been saved.

7. A screen appears that lets you review your basal pattern. Press **✓** to review all the settings. If you need to make any changes, press **⬅** to return to the edit screen.

**Note:** If you press **⬅** to return to the edit screen without saving, your changes will not be saved. After making any changes on the edit screen, select **Review** to return to the review screen.

8. Select **Save**. If you do not select **Save**, your changes are not saved.



**CAUTION:** If you have not selected **Save** after settings are entered and the screen goes dark, the entered settings will not be saved.










**Turn On and Set Up the Bolus Wizard Feature**

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

## Menu > ⚙ > Delivery Settings > Bolus Wizard Setup

1. Select **Bolus Wizard** to turn the feature on.
2. Press **✓** to see the rest of the message then select **Next**.
3. Review the description of carb ratio and select **Next**.
4. Press **○** to enter **End** time. Press **○** again to leave the end time at **12:00A**. If you are setting multiple carb ratios use **^** or **✓** to change the end time, then press **○**.
5. Press **^** or **✓** to enter **g/U** then press **○**. Continue adding end times and carb ratios as needed.
6. Select **Next** when done.
7. Review the description of insulin sensitivity factor and select **Next**.
8. Press **○** to enter **End** time. Press **○** again to leave the end time at **12:00A**. If you are setting multiple insulin sensitivity factors, use **^** or **✓** to change the end time, then press **○**.
9. Press **^** or **✓** to enter **mg/dL per U**, then press **○**. Continue adding end times and insulin sensitivity factors as needed.
10. Select **Next** when done.
11. Review the description of BG target and select **Next**.
12. Press **○** to enter **End** time. Press **○** again to leave the end time at **12:00A**. If you are setting multiple BG targets, use **^** or **✓** to change the end time, then press **○**.
13. Press **^** or **✓** to enter the low limit (Lo) and then press **○**.
14. Press **^** or **✓** to enter the high limit (Hi) and then press **○**. Continue adding end times and BG targets as needed.
15. Select **Next** when done.
16. Review the description of active insulin time and select **Next**.
17. Press **○** to enter the duration for the active insulin time.
18. Press **^** or **✓** to enter the hours of duration then press **○**.
19. Select **Save**.

The Bolus Wizard feature setup is now complete.





## Turn Sensor On

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

## Menu > ⚙ > Device Settings

- Select **Sensor** to turn the feature on. You can now access the **SmartGuard** menu and enter the SmartGuard settings.





## Pair the Transmitter

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.


### Menu >

1. Attach Guardian™ 4 transmitter to the charger and place the transmitter (attached to the charger) next to the pump.
2. Select **Pair New Device** on the pump and immediately remove the transmitter from the charger.
3. The **Select Device** screen appears with a list of available devices. Select the CGM device that matches the serial number indicated on the back of the transmitter.

Your pump displays a message if the pump and transmitter are paired successfully. When the transmitter is communicating with the pump, the Sensor feature is turned on and  appears on the Home screen.
















SN XXXXXXXXXXXX




## Programming High SG Settings

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

### Menu > > Alert Settings

1. Select **High Alert**. The **High Setup** screen appears.
2. Press  on the time segment. Press  again to leave the **End** time at **12:00A**. If setting multiple time segments, press  to enter the end time of the first segment, and press .
3. Press  or  to set the high limit (Hi) and press . Press  to continue.
4. Select each setting you wish to turn on. If a setting is on, select it again to turn it off.
5. Once settings are selected, select **Next**. If needed, continue adding time segments and settings to complete the 24-hour period. When done, select **Review**.
6. Verify that settings are correct and select **Save**.
7. If snooze time needs to be changed, press  to **Snooze High & Low** and press . Select **High Snooze**. Press  or  to set up the desired duration and press . Select **Save**.

Your High Setup is now complete.







## Programming Low SG Settings

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

## Menu > ⚙ > Alert Settings

1. Select **Low Alert**. The **Low Setup** screen appears.
2. Press **○** on the time segment. Press **○** again to leave the **End** time at **12:00A**. If setting multiple time segments, press **^** to enter the end time of the first segment, and press **○**.
3. Press **^** or **v** to set the low limit (Lo) and press **○**. Press **○** to continue.
4. Select each setting you wish to turn on. If a setting is on, select it again to turn it off.
5. Once settings are selected, select **Next**. If needed, continue adding time segments and settings to complete the 24-hour period. When done, select **Review**.
6. Verify that settings are correct and select **Save**.
7. If snooze time needs to be changed, press **v** to **Snooze High & Low** and press **○**. Select **Low Snooze**. Press **^** or **v** to set up the desired duration and press **○**. Select **Save**.

Your Low Setup is now complete.




## Turn SmartGuard On

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

## Menu > ☺

1. Select **SmartGuard** to turn the SmartGuard feature on.
2. Review the message on the screen and select **OK**.
3. Review the SmartGuard Checklist to see if any additional requirements must be met before the SmartGuard feature can be activated. Press **v** to view the checklist.







**Note:** The SmartGuard warm-up period begins at the first midnight after the pump starts delivering insulin. If you turn SmartGuard on before the warm up is complete, you may receive an alert at midnight to enter a BG for SmartGuard. You may want to wait until day 3 of your pump use to avoid this alert at midnight.


## Programming SmartGuard Settings

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

**Note:** Make sure you consult with your healthcare provider on your personalized settings before programming the SmartGuard settings.

Menu >  > SmartGuard Settings

1. Review the **Target**.
2. Press  to change the target.
3. Press  or  to the desired target and then press .
4. Check to make sure that **Auto Correction** is set to **On**. Select **Save**.



**Note:** Please take some time to review ALL your settings and ensure they have been correctly saved.

**Technical Support**

**In Need of Technical Assistance?**

24-Hour Technical Support is here to help

Call 1.800.646.4633

For medical advice, contact your healthcare team.



# Medtronic



**Medtronic MiniMed**

18000 Devonshire Street  
Northridge, CA 91325

USA

1 800 646 4633

+1 818 576 5555

[www.medtronicdiabetes.com](http://www.medtronicdiabetes.com)

RF: M994838A001

© 2023 Medtronic  
M037928C001\_1  
2023-05-10



\*M037928C001\*