Français

Mode d'emploi

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


English

Indications for Use

The Sof-set is indicated for the subcutaneous infusion of medication, including insulin, from a Paradigm infusion pump and reservoir.

Contraindications

The Sof-set family of infusion sets is indicated for subcutaneous use only and not for IV (intravenous) infusion or the infusion of blood or blood products.

Warning

Inaccurate medication deliver, infection and/ or site irritation may result from improper insertion and maintenance of infusion site. Remove needle guard before inserting infusion set.

If using this infusion set for the first time, do first set-up in the presence of healthcare provider.

Since different individuals have varying amounts of subcutaneous fat, be sure to choose appropriate cannula length to fit your needs.

Do not leave air in infusion set. Prime completely.

Check frequently to make sure soft cannula remains firmly in place as you may not feel pain if it pulls out. The soft cannula must always be completely inserted to receive full amount of medication.

If infusion site becomes inflamed, replace set, and use new site until first site has healed.

Replace the Sof-set if the tape becomes loose, or if the soft cannula becomes fully or partially dislodged from the skin.

Replace infusion set every 48-72 hours according to Center for Disease Control guidelines, or per healthcare provider's instructions.

Do not reinsert introducer needle into infusion set. Re-insertion could cause tearing of soft cannula which would result in unpredictable medication flow.

Never prime set or attempt to free a clogged line while set is inserted. You may accidentally inject too much medication.

Do not put disinfectants, perfumes, or deodorants on infusion set as these may affect integrity of set.

Dispose of infusion set and needle safely, in sharps container, after single use. Do not clean or re-sterilize.

Infusion Set • Dispositif de perfusion • Infusionsset • Equipo de infusión • Infusieset • Set di infusione • Slangeset • Infusioipakkkaus • Infusionsset • Infusionssæt • Infúzní set • Zestaw infuzyjny • Σετ έγχυσης

^[1] Infusion Set • Dispositif de perfusion • Infusionsset • Equipo de infusión • Infusieset • Set di infusione • Slangeset • Infusioipakkkaus • Infusionsset • Infusionssæt • Infúzní set • Zestaw infuzyjny • Σετ έγχυσης

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