


MiniMed™ 780G system Detailed Menu Map



<p>1 </p> <p>Insulin</p> <p>Bolus/Stop Bolus*</p> <ul style="list-style-type: none"> • Bolus Wizard • Manual Bolus • Preset Bolus* • Delivery Settings ➔ <p>Basal/Cancel Temp Basal*</p> <ul style="list-style-type: none"> • Temp Basal • Preset Temp* • Basal Patterns • Basal Pattern Setup ➔ <p>Suspend All Delivery/Resume Basal*</p> <p>Delivery Settings ➔</p>	<p>2 </p> <p>History & Graph</p> <p>History</p> <ul style="list-style-type: none"> • Summary • Daily History • Alarm History • Paired Sensors <p>Sensor Glucose Review</p> <ul style="list-style-type: none"> • High Limit • Low Limit • Days to Average <p>Graph</p> <p>Time in Range</p>	<p>3 </p> <p>SmartGuard</p> <p>SmartGuard Checklist</p> <p>Temp Target/Cancel Temp Target*</p> <p>SmartGuard Settings</p> <ul style="list-style-type: none"> • Target • Auto Correction <p>SmartGuard (On and Off)</p>
<p>6 </p> <p>Settings</p> <p>Alert Settings</p> <ul style="list-style-type: none"> • High Alert • Low Alert • Snooze High & Low • Reminders <p>Delivery Settings</p> <ul style="list-style-type: none"> • Bolus Wizard Setup • Basal Pattern Setup • Max Basal/Bolus • Dual/Square Wave • Bolus Increment • Bolus Speed • Preset Bolus Setup • Preset Temp Setup 	<p>6 </p> <p>Settings continued</p> <p>Device Settings</p> <ul style="list-style-type: none"> • Sensor (On and Off) • Time & Date • Device Info • Display • Block Mode • Self Test • Review Settings • Manage Settings • Easy Bolus • Auto Suspend • Language 	<p>4 </p> <p>Sound & Vibration</p> <p>Silence Sensor Alerts</p> <p>Volume</p> <p>Sound</p> <p>Vibration</p> <p>Alert Settings ➔</p> <p>5 </p> <p>Blood Glucose</p> <p>BG</p> <p>7 </p> <p>Status</p> <p>Stop Bolus</p> <p>Suspend All Delivery/Resume Basal*</p> <p>SmartGuard Checklist</p> <p>Pump</p> <p>Sensor</p> <p>8 </p> <p>Paired Devices</p> <p>Pair New Device</p> <p>Pair CareLink</p> <p>Mobile</p> <p>Meter</p> <p>CGM</p> <p>9 </p> <p>Reservoir & Set</p> <p>New Reservoir & Set</p> <p>New Reservoir Only</p> <p>New Set Only</p> <p>Fill Cannula</p>